

MODIFICATION OF POLYMER FOILS WITH EUV (EXTREME-ULTRAVIOLET) RADIATION FOR APPLICATIONS IN BIOMEDICAL TECHNOLOGY

Project partners:

Institute of Applied Physics, Johannes-Kepler University Linz, Linz, Austria (Coordinator)
 Institute of Optoelectronics, Military University of Technology, Warsaw, Poland
 Rigaku Innovative Technology Europe, Prague, Czech Republic
 Institute of Biophysics, Johannes Kepler University Linz, Linz, Austria



Scientific and Technological Project Goals:





Novel methods and techniques of polymer surface modification for biomedical applications using extreme ultraviolet (EUV) radiation has been developed and a new compact laser plasma EUV source based on a gas puff target dedicated for this technique was built.

LASER PLASMA EUV SOURCE BASED ON A GAS PUFF TARGET

Laser plasma EUV source

Double-stream gas puff target



Laser plasma EUV source for modification polymer surfaces



BIOCOMPATIBILITY STUDIES

Cultivation of CHO cells at EUV irradiated PET

SEM



Phase Contrast Microscope



- good alignment of Chinese hamster ovary (CHO) cells along the direction of the walls,
- CHO cells showed only bad adhesion at the irradiated surfaces and no alignment for samples irradiated with UV laser at 193 nm,
- chemical surface modification is more pronounced for EUV irradiation.

PUBLICATIONS

- A. Bartnik, H. Fiedorowicz, R. Jarocki, J. Kostecki, M. Szczurek, P.W. Wachulak, „Laser-plasma EUV source dedicated for surface processing of polymers” Nucl. Instr. Meth. Phys. Res. A 647 (2011) 125
- A. Bartnik, H. Fiedorowicz, R. Jarocki, J. Kostecki, M. Szczurek, O. Chernyayeva, J.W. Sobczak, „EUV-induced physico-chemical changes in near-surface layers of polymers” J. Electron Spectrosc. Relat. Phenom. 184 (2011) 270
- A. Bartnik, H. Fiedorowicz, S. Burdyńska, R. Jarocki, J. Kostecki, M. Szczurek, „Combined effect of EUV irradiation and acetone treatment on PET surface”, Appl. Phys. A 103 (2011) 173
- B. Reisinger, M. Fahrner, I. Frischauf, S. Yakunin, V. Svorcik, P. Sajdl, H. Fiedorowicz, A. Bartnik, C. Romanin, J. Heitz, „EUV micro-patterning for biocompatibility control of PET”, Appl. Phys. A 100 (2010) 511
- A. Bartnik, H. Fiedorowicz, R. Jarocki, J. Kostecki, M. Szczurek, A. Biliński, O. Chernyayeva, J.W. Sobczak, „Physical and chemical modifications of PET surface using a laser-plasma EUV source”, Appl. Phys. A 99 (2010) 831
- A. Bartnik, H. Fiedorowicz, R. Jarocki, J. Kostecki, M. Szczurek, „PMMA and FEP surface modifications induced with EUV pulses in two selected wavelength ranges”, Appl. Phys. A 98 (2010) 61
- A. Bartnik, H. Fiedorowicz, R. Jarocki, J. Kostecki, R. Rakowski, A. Szczurek, M. Szczurek, „Micro- and Nanoprocessing of Polymers Using a Laser Plasma Extreme Ultraviolet Source”, Acta Physica Polonica A, 117 (2010) 384
- A. Bartnik, H. Fiedorowicz, R. Jarocki, J. Kostecki, A. Szczurek, M. Szczurek, „Ablation and surface modifications of PMMA using a laser-plasma EUV source” Appl. Phys. B 96 (2009) 727
- A. Bartnik, H. Fiedorowicz, R. Jarocki, J. Kostecki, A. Szczurek, M. Szczurek, „Creation of nanostructures on polymer surfaces irradiated with extreme ultraviolet pulse” Acta Phys. Pol. A 116 (2009) S108

MODIFICATION OF POLYMER SURFACES USING EXTREME ULTRAVIOLET (EUV)

Typical SEM images of EUV irradiated polymers

